

Information for opening GUI:

1. in terminal type: `python particle_label_gui_v9.py`
2. A window showing files on your computer should pop up.
3. Within the file window navigate to the subfolders of the “**particle_recognition**” folder, each subfolder should contain a folder called “**adjustedPNG2**”
4. Click “**adjustedPNG2**” , you should see a list one png files, click open button in bottom right corner.
5. The first image to label should now be displayed

NOTE: GUI will always start with first image in a file. If images have already been labeled but GUI is opened on the same file a second time the labels will be overwritten! Once you open a file to label make sure to finish all images before quitting GUI.

Information for labeling images:

- Atomic columns must be resolved to either label particle “yes” or “no”
- If atomic columns are not resolved particle must be labeled “null”
- Particles containing stacking fault should be labeled “yes”
- Particles with atomic columns resolved but no stacking fault should be labeled “no”
- Circle should be centered on particle and as tight fitting as possible without cutting out any part of the particle

Key codes for running GUI:

- “*SHIFT + CLICK*” = select particle
- “*ARROW KEYS*” = move circle to center circle on particle
- “*e*” = enlarge size of circle
- “*d*” = decrease size of circle
- “*y*” = label particle YES
- “*n*” = label particle NO
- “*o*” = label particle NULL
- “*ENTER*” = save labels for this image and load next image, once all files have been labeled hitting enter will quit GUI
- “*q*” = quit GUI